



MEDELEQ FRESCO SILICONE LIST 2005

IMAGE	PRODUCT CODE	PRODUCT DESCRIPTION	SIZE	QTY	
PADDING AND PROTECTION PRODUCTS					
	F-00001-01	Tubifoam	Small	Each	
	F-00001-02	Tubifoam	Medium	Each	
	F-00001-03	Tubifoam	Large	Each	
	F-00001-04	Tubifoam	X-Large	Each	
	F-00003-01	Foam Toe Separators	Small	4 Per Bag	
	F-00003-02	Foam Toe Separators	Medium	4 Per Bag	
	F-00003-03	Foam Toe Separators	Large	4 Per Bag	
	F-00003-04	Foam Toe Separators	X-Large	3 Per Bag	
	F-00004-01	Toe Prop Leather	Small	Pair	
	F-00004-02	Toe Prop Leather	Medium	Pair	
	F-00004-03	Toe Prop Leather	Large	Pair	
	F-00006-00	Corn Protector - Latex Foam	One Size	Each	
	F-00014-00	Heel Raiser (Foam) with Spur Cut out – 6mm	XS, S, M	Pair	
	F-00018-01	Method Carré - Treatment pack for ingrowing toenails		Box of 20 units	
	F-00020-01	Halux Valgus Night Splint	Small	Pair	
	F-00020-02	Halux Valgus Night Splint	Medium	Pair	
	F-00020-03	Halux Valgus Night Splint	Large	Pair	
	F-00021-01	Moleskin Toe Separators	Small	4	
	F-00021-02	Moleskin Toe Separators	Medium	4	
	F-00021-03	Moleskin Toe Separators	Large	4	
	F-00021-04	Moleskin Toe Separators	X-Large	3	










IMAGE

PRODUCT CODE

PRODUCT DESCRIPTION

SIZE

QTY

		FRESIL SILICONE INSOLE RANGE			
	F-00031-01	Comfort Heel Raisers – one density – 4mm	Small	Pair	
	F-00031-02	Comfort Heel Raisers – one density	Medium	Pair	
	F-00031-03	Comfort Heel Raisers – one density	Large	Pair	
	F-00031-04	Comfort Heel Raisers - two densities with softer blue dot – 4mm	Small	Pair	
	F-00031-05	Comfort Heel Raisers - two densities with softer blue dot	Medium	Pair	
	F-00031-06	Comfort Heel Raisers - two densities with softer blue dot	Large	Pair	
	F-00037-01	Heel Cup - two densities with softer blue dot	Small - 9.1cm	Pair	
	F-00037-02	Heel Cup - two densities with softer blue dot	Medium - 10.9cm	Pair	
	F-00037-03	Heel Cup - two densities with softer blue dot	Large - 12cm	Pair	
	F-00038-01	3/4 length silicone insoles with raised arch & softer blue dot in heel	Small - 16.4cm	Pair	
	F-00038-02	3/4 length silicone insoles with raised arch & softer blue dot in heel	Medium - 17cm	Pair	
	F-00038-03	3/4 length silicone insoles with raised arch & softer blue dot in heel	Large - 17.7cm	Pair	
	F-00040-05	Full Length silicone insoles with two densities with metatarsal elevation	X/Small - Length 23.5cm	Pair	
	F-00040-06	Full Length silicone insoles with two densities with metatarsal elevation	Small - Length 25cm	Pair	
	F-00040-07	Full Length silicone insoles with two densities with metatarsal elevation	Medium - Length 26.3cm	Pair	
	F-00040-08	Full Length silicone insoles with two densities with metatarsal elevation	Large - Length 27.7cm	Pair	
	F-00040-09	Full Length silicone insoles with two densities with metatarsal elevation	X/Large - Length 29cm	Pair	
	F-00071-00	Silicone met protector	One size	Pair	
					

IMAGE

PRODUCT CODE

PRODUCT DESCRIPTION

SIZE

QTY



FELTS AND MOLESKINS



F-00506-01	Moleskin Fresco colour beige 100% Cotton - hypoallergenic adhesive - 0.5mm thick	Roll 25 cm x 3m		
F-00506-02	Moleskin Fresco colour pink 100% Cotton hypoallergenic adhesive- 0.5mm thick	Roll 25 cm x 3m		
F-00507-01	Felt 1mm 50% wool and 50% rayon viscose with hypoallergenic adhesive	15 sheets of 30 x 20 cm	One Box	
F-00507-02	Felt 2mm 50% wool and 50% rayon viscose with hypoallergenic adhesive	10 sheets of 30 x 20 cm	One Box	
F-00507-03	Felt 3mm 50% wool and 50% rayon viscose with hypoallergenic adhesive	6 sheets of 30 x 20 cm	One Box	
F-00507-04	Felt 4mm 50% wool and 50% rayon viscose with hypoallergenic adhesive	6 sheets of 30 x 20 cm	One Box	
F-00507-06	Felt 6mm 50% wool and 50% rayon viscose with hypoallergenic adhesive	4 sheets of 30 x 20 cm	One Box	
F-00507-08	Felt 8mm 50% wool and 50% rayon viscose with hypoallergenic adhesive	3 sheets of 30 x 20 cm	One Box	
F-00508-00	Velvet Moleskin Rubber adhesive hypoallergenic 100% cotton 0.8mm thick	15 sheets of 30 x 20 cm	One Box	
F-00508-10	Velvet Moleskin Rubber adhesive hypoallergenic 100% cotton 0.8mm thick	32cm Wide	Price Per Metre	



F-00509-00	Moleskin Fresco Sport 100% cotton hypoallergenic adhesive approx 0.8mm thick	15 sheets of 30 x 20 cm	One Box	
F-00509-10	Moleskin Fresco Sport 100% cotton hypoallergenic adhesive approx 0.8mm thick	32cm Wide	Price Per Metre	
F-00514-00	Moleskin Fresco Sport 100% cotton hypoallergenic adhesive approx 2mm thick	10 sheets of 30 x 20 cm	One Box	
F-00607-01	Fresco Prote – White Foam with Hypoallergenic Adhesive 30 x 20cm per box	4 Sheets	One Bag	



IMAGE

PRODUCT CODE

PRODUCT DESCRIPTION

SIZE

QTY



MISCELLANEOUS PRODUCTS



	F-00910-00	Martin goniometric ruler For measuring different angles of the foot. It has a 170mm reglet. Measures: 200 x 110mm	One	One	
	F-00920-00	Goniometric reglet Of transparent plastic and can be folded up. It has graduated segments, Concentric circles, graduated semi-circle, inverted T, central graduated circle and lateral marks. Measures length: 360mm, closed 205mm	One	One	
	F-00921-00	Goniometric tablet Perthes Type. For heel measures. Measures: 102 x 120mm	One	One	
	F-01111-00	Podotrack (100 sheets) Foot Impression Sheets No need for tints, liquids etc. Backing sheet allows annotation	100 sheets	One Box	
	F-01551-00	Podoscope Acrylic with Metal, Fluorescent Light	One	One	
	F-00059-00	Skeleton of the foot for practicing procedures	One	One	
	F-00065-00	Poster - Artists sketch drawing of the foot in 4 positions Coloured Width: 100cm, High: 35.5cm	111.4 x 35.5cm	One	
	F-00080-00	Bronze figure of boy examining his foot Measures: height 195mm, base: 150mm, weight: 2800gr	Weight 2,800 grams	One	

IMAGE

PRODUCT
CODE

PRODUCT DESCRIPTION

SIZE

QTY



SILICONES AND CATALYSTS



F-01903-00	Spatula Inox - for liquid silicones' working		One	
F-01904-00	Rigid work surface board - for liquid silicones'.	20 x 25 cm.	One	
F-01906-01	Polymer of silicone - Liquid silicone - semi-rigid - medium elasticity Used to repair silicone orthotics, or as a joining material with other silicones.	250 gram	One	
F-01907-01	Silicone RTV-11504 - Liquid silicone - Very elastic - flexible - with a high adherence	250 gram	One	
F-01908-00	Reaktol - liquid catalyst (for use with Bland Rose and Fresco Silicones)	20ml	One	
F-01910-00	Fresco catalyst - Red paste for mixing all types of Fresco Silicones. Its red colour gives a soft red colour to the final orthosis. It is a little bit slower and gives a softer result orthosis than the liquid catalyst.	20 gram	One	
F-01920-00	Silicone oil Used to mix with all the Fresco Silicones	200 gram	One	
F-01921-01	Silicone 1.400 - Putty - grey colour Hard Silicone, semi rigid with grey colour. Semi-rigid silicone elastomer in the form of a light grey coloured paste. Quick setting	500 gram	One	
F-01922-02	Silicone Orthesil Putty - orange colour. Used for rigid orthoses. This silicone is not elastic but has a long life span. It is recommended to be used for the correction or re-alignment of toes.	500 gram	One	

IMAGE

PRODUCT CODE

PRODUCT DESCRIPTION







SIZE

QTY



SILICONES AND CATALYSTS



	F-01923-00	Silicone Fresco - Putty - semi-soft silicone, soft to the touch, with average elasticity.	500 gram	One	
	F-01924-00	Silicone DUO 10 – A + B silicone – Hardness Shore A 10	250 gr + 250 gr.	One pack (500 gr)	
	F-01926-00	Starting Kit – consists of: 100 gr. Bland-Rosé + 100 gr. Fresco + 1 Reaktol (Catalyst)		One	
	F-01930-01	Silicone Bland-Rose Includes 1 x Reaktol Catalyst Soft silicone with high elasticity. The ideal silicone for palliative orthosis.	500 gram	One	
	F-01933-00	Silicone Bland-Nature - Putty - Pink skin colour Used as a moulding material that produces vulcanised silicone of minimum hardness with excellent mechanical resistance.	500 gram	One	
	F-01940-01	Sorbopres - Silicone sheet with adhesive back - Medium hardness: 15° Shore	2mm thick 30 x 27.5cm	One	
	F-01940-02	Sorbopres - Silicone sheet with adhesive back - Medium hardness: 15° Shore	3mm thick 30 x 27.5cm	One	
	F-01941-01	Fresco Pad-A - Silicone film - semi adhesive- long durability Shore A1	1mm thick 22 x 12cm	One	
	F-01941-04	Fresco Pad - Silicone sheet - semi adhesive Silicone Sheet used as a protector for any part of the foot. The fabric in the middle gives a longer durability to the product and keeps the hardness very low.	2mm thick 30 x 11cm	One	

IMAGE

PRODUCT CODE

PRODUCT DESCRIPTION









SIZE

QTY



FRESCO SILICONE PRODUCT RANGE



	F-00033-01	Pure gel cap	1.5cm Diameter S/M	Each	
	F-00033-03	Pure gel cap	2.0cm Diameter L/XL	Each	
	F-00033-02	Digital Cap Ribbed Fabric	1.5cm Diameter S/M	Each	
	F-00033-04	Digital Cap Ribbed Fabric	2cm Diameter L/XL	Each	
	F-00036-01	Gel Sheet Covered one-side with cloth	1mm thick 15 x 20cm	1 Sheet	
	F-00036-02	Gel Sheet Covered one-side with cloth	2mm thick 15 x 20cm	1 Sheet	
	F-00036-03	Gel Sheet Covered both sides	2mm thick 15 x 20cm	1 Sheet	
	F-00041-01	Pure gel ring	1.5cm Diameter S/M	Each	
	F-00041-02	Pure gel ring	2cm Diameter L/XL	Each	
	F-00041-04	Digital Ring Pad/Corn Pad	Mini	Each	
	F-00041-05	Digital Ring Pad/Corn Pad	Small/Medium	Each	
	F-00041-06	Digital Ring Pad/Corn Pad	Large/X-Large	Each	
	F-00042-01	Hallux Valgus Gel Protector (Fabric Housing) Gel on H/V Joint only	Small/Medium	Pair	
	F-00042-02	Hallux Valgus Gel Protector (Fabric Housing) Gel on H/V Joint only	Large/X-Large	Pair	
	F-00042-03	Hallux Valgus Gel Protector (Fabric Housing) Gel coated throughout	Small/Medium	Pair	
	F-00042-04	Hallux Valgus Gel Protector (Fabric Housing) Gel coated throughout	Large/X-Large	Pair	
	F-00042-07	Small Bunion Gel Protector	One size	Each	
	F-00043-01	Metatarsal Gel Strap	8cm width (Elastic will stretch) Small/Medium	Pair	
	F-00043-02	Metatarsal Gel Strap	10cm width (Elastic will stretch) Large/X-Large	Pair	

IMAGE

PRODUCT CODE

PRODUCT DESCRIPTION

SIZE

QTY




FRESCO SILICONE PRODUCT RANGE



F-00044-01	Metatarsal Protector	Small/Medium	Pair	
------------	----------------------	--------------	------	--

F-00044-02	Metatarsal Protector	Large/X-Large	Pair	
------------	----------------------	---------------	------	--



F-00045-03	Digital Gel Dot Flat on adhesive sheet		9 Dots	
------------	--	--	--------	--

F-00045-04	Digital Gel Dots Concave on adhesive sheet		9 Dots	
------------	--	--	--------	--



F-00046-00	Pure Gel Square – No back	10cm x 10cm	2 Sheets	
------------	---------------------------	-------------	----------	--



F-00047-00	Pure Gel Square – With adhesive back	10cm x 10cm	2 Sheets	
------------	--------------------------------------	-------------	----------	--



F-00048-01	Digistrip - Roll	72cm x 15mm Small	1 Roll	
------------	------------------	-------------------	--------	--

F-00048-02	Digistrip - Roll	72cm x 20mm Medium	1 Roll	
------------	------------------	--------------------	--------	--

F-00048-03	Digistrip - Roll	72cm x 25mm Large	1 Roll	
------------	------------------	-------------------	--------	--



F-00050-01	Gel Tube – Mesh Material	1.5cm x 15cm S/M	Each	
------------	--------------------------	------------------	------	--

F-00050-02	Gel Tube – Mesh Material	2cm x 15cm L/XL	Each	
------------	--------------------------	-----------------	------	--

F-00050-03	Gel Tube – Ribbed Material	1.5cm x 15cm S/M	Each	
------------	----------------------------	------------------	------	--



F-00050-04	Gel Tube – Ribbed Material	2cm x 15cm Large	Each	
------------	----------------------------	------------------	------	--

F-00050-05	Gel Tube – Ribbed Material	2.5cm x 15cm X/Large	Each	
------------	----------------------------	----------------------	------	--



F-00051-01	Pure gel tube	1.5 cm x 8 cm S/M	Each	
------------	---------------	-------------------	------	--

F-00051-02	Pure gel tube	2 cm x 8 cm L/XL	Each	
------------	---------------	------------------	------	--

IMAGE

PRODUCT CODE

PRODUCT DESCRIPTION




SIZE

QTY



FRESCO SILICONE PRODUCT RANGE



	F-00052-01	Gel Half Sock (gel coating on sole: 3mm thick) Ideal for geriatrics, diabetics and feet with no plantar fat	Small	Pair	
	F-00052-02	Gel Half Sock (gel coating on sole)	Medium	Pair	
	F-00052-03	Gel Half Sock (gel coating on sole)	Large	Pair	
	F-00052-04	Gel Half Sock (gel coating on sole)	X/Large	Pair	
	F-00053-01	Gel Insole – (Fabric coated both sides: 2mm thick gel film)	Small	Pair	
	F-00053-02	Gel Insole – (Fabric coated both sides)	Medium	Pair	
	F-00053-03	Gel Insole – (Fabric coated both sides)	Large	Pair	
	F-00054-01	Toe Spreader (Shape – Thick)	Small/Med	Each	
	F-00054-02	Toe Spreader (Shape – Thick)	Large/XL	Each	
	F-00055-01	Gel Separator (Shape - Shark Fin)	Mini	Each	
	F-00055-02	Gel Separator (Shape - Shark Fin)	Small/Med	Each	
	F-00055-03	Gel Separator (Shape - Shark Fin)	Large	Each	
	F-00056-01	Pure Gel Toe Props	Small/Med	Pair	
	F-00056-02	Pure Gel Toe Props	Large/XL	Pair	
	F-00057-01	Pure Gel Discs – one size only		Each	
	F-00058-01	Gel Foot Cover	Small/Med	Pair	
	F-00058-02	Gel Foot Cover	Large/XL	Pair	
	F-00042-05	Hallux Valgus Protector – Pure Gel Bunion Shield	One Size	Each	



GEL PROPERTIES

All Fresco gel products contain a cross-linked three-dimensional polymer network swollen in a medical grade mineral oil. The gel contains a medical grade mineral oil that softens, moisturises and lubricates the skin. The oil is slowly released to break down keratotic lesions and reduce and prevent keloid and hypertrophic scarring.

- Friction to the skin is greatly reduced because the top and bottom surfaces of the gel move independently of each other.
- Because the gel flows in all directions when subjected to weight or pressure, shock, vibration and shear forces are minimised.
- Unlike standard felts and foams the polymer gel will not bottom out.
- The gel conforms to bony prominences giving superb comfort and pain relief.
- The gel has a “memory” and will quickly return to its original form after being stretched or manipulated. It has an elasticity of 400%.
- All products in the range are washable and reusable unlike many other footcare products on the market today.
- Washable and reusable.
- Fresco polymer gel is hypoallergenic, non-toxic and does not support bacterial growth.



Distributed By:
MEDELEQ PTY LTD
Ph: 07 55 962422
E: sales@medeleq.com.au
W: www.medeleq.com.au

