

















Medeleq Baby Scales Comparison

Model	M118600-01	M118600	M112600	MED-BAK	MED-BARL	M108800	M111600
							
Type	Precision Digital	Precision Digital	Digital	Digital	Digital	Mechanical	Digital - Hand Held
Capacity	20KG	20KG	20KG	20KG	20KG	20KG	20KG
Graduations	5g	5g	5g	5g	5g	50 g < 10 kg > 100 g	20g
Display	Backlit LCD	Backlit LCD	Backlit LCD	Backlit LCD	Backlit LCD	Dial Gauge with 2 needles	Backlit LCD
Power Source	4 x AA Battery	4 x AA Battery	4 x AA Battery	4 x AA Battery	4 x AA Battery	Not Required	2 x CR2032
Length Measurement	Yes- Digital Calipers	No	No	No	Yes - Mechanical Calipers	No	No
Length Range	400 - 800mm, 2mm Grad				350 - 700mm, 1mm Grad		
Dimensions	644 x 350 x 135 mm	600 x 350 x 135 mm	600 x 385 x 85 mm	580 x 325 x 120 mm	600 x 400 x 120 mm	550 x 259 x 343 mm	
Weight	3.4 kg	3.2 kg	2.6 kg	2.3 kg	2.5 kg	4.5 kg	200g
Features	The optionally integrated digital length measuring device will take precise readings of infants up to 80 cm in 2 mm steps. The infant can be freely positioned inside the weighing tray, no additional handling necessary. The measuring results will be shown on the LCD display of the baby weighing scale. All weighing functions and specifications are equivalent to ADE M118600.	Hospital Baby Scale. The gently contoured baby tray offers generous as well as secure space for small infants and even for seated toddlers. The fine graduation of 5 g up to 20 kg, together with the auto-hold function and an electronic dampening system guarantees precise and fast weighing results. The clear backlit display ensures easy reading also in dimmed surroundings. Simple to operate by five clearly marked soft touch buttons.	Large, but lightweight baby weighing scale with a generous weighing tray, offering plenty of space to infants and toddlers. Ideally suited for GPs, clinics and for home visiting nurses as a portable scale. A clear 30 mm display and the simple three key operation make the scale perfectly easy to use.	Baby / Toddler weighing scale. LCD display with backlight. Weight saving function. Last weight recall function. Kg/lb units exchange. 4 x AA batteries power. Weighing Tray can be removed to enable toddlers to stand on the device.	Shallow Tray Baby weighing scale. LCD display with backlight. Weight saving function. Last weight recall function. Kg/lb units exchange. 4 x AA batteries power. Integrated length measuring Rod with calipers.	Quick and easy to use baby and toddler scale. Clear indication of the weight by two red revolving indicator needles. Robust metal housing. For general use and also for field workers in rural situations.	Electronic hanging type baby scale – small and practical device made of steel, with large hand grip and hook for weighing trouser. Weighs babies up to 20 kg with a precision of 20 g. Two chrome plated hooks are included.
	Auto Zero Setting, Tare, Hold, Low Batt, Auto Off, Weight or Length Display High End Load Cells	Auto Zero Setting, Tare, Hold, Low Batt, Auto Off High End Load Cells	Auto Zero Setting, Tare, Hold, Low Batt, Auto Off	Auto Zero Setting, Tare, Hold, Low Batt, Auto Off, Last Weight recall.	Auto Zero Setting, Tare, Hold, Low Batt, Auto Off, Last Weight recall.	Manual zero setting, Tare	Auto Zero Setting, Tare, Hold, Low Batt
Accessories	 Carry Bag MZ10060	 Carry Bag MZ10060	 Carry Bag MZ10060	 Carry Bag MZ10060	 Carry Bag MZ10060	 Carry Bag MZ10060	 Weighing Sling MZ111600-002
	 Mains Adaptor	 Mains Adaptor					